

PE anti-human CD275 (B7-H2, ICOSL)

Mouse IgG2b, kappa
2D3



BioLegend
From Revvity

25 tests

Artikel-Nr.: B309403 | **BioLegend** | Hersteller-Nr.: 309403

94,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 2D3
Format: PE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE anti-human CD275 (B7-H2, ICOSL)

Clone: 2D3
Isotype: Mouse IgG2b, kappa
Other_Names: CD275, ICOSL, B7H2, ICOS-L, B7RP-1, B7-RP1
Size: 25 tests

Reactivity: Human
Format: PE
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1515>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>