

Purified anti-human CD282 (TLR2)

Mouse IgG2a, kappa
TL2.1



BioLegend
From Revvity

100 µg
Artikel-Nr.: B309702 | **BioLegend** | Hersteller-Nr.: 309702

207,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: TL2.1
Format: PURE
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD282 (TLR2)

Clone: TL2.1
Isotype: Mouse IgG2a, kappa
Other_Names: Toll like receptor 2, TIL4, CD282
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1502>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>