

## Alexa Fluor® 488 anti-human CD282 (TLR2)

Mouse IgG2a, kappa  
TL2.1



**BioLegend**  
From Revvity

100 tests

Artikel-Nr.: B309712 | BioLegend | Hersteller-Nr.: 309712

**235,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FC, ICFC  
**Clone:** TL2.1  
**Format:** A488  
**Isotype:** Mouse IgG2a, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 Tests

Alexa Fluor® 488 anti-human CD282 (TLR2)

Clone: TL2.1

Isotype: Mouse IgG2a, kappa

Other\_Names: Toll like receptor 2, TIL4, CD282

Size: 100 tests

Reactivity: Human

Format: A488

App\_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3148>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>