

Alexa Fluor® 647 anti-human CD282 (TLR2)

Mouse IgG2a, kappa
TL2.1



BioLegend
From Revvity

100 tests

Artikel-Nr.: B309714 | BioLegend | Hersteller-Nr.: 309714

235,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: TL2.1
Format: A647
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

Alexa Fluor® 647 anti-human CD282 (TLR2)

Clone: TL2.1

Isotype: Mouse IgG2a, kappa

Other_Names: Toll like receptor 2, TIL4, CD282

Size: 100 tests

Reactivity: Human

Format: A647

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3149>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>