

Purified anti-human CD137 (4-1BB)

Mouse IgG1, kappa
4B4-1



BioLegend
From Revvity

100 µg

Artikel-Nr.: B309802 | **BioLegend** | Hersteller-Nr.: 309802

105,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 4B4-1
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD137 (4-1BB)

Clone: 4B4-1

Isotype: Mouse IgG1, kappa

Other_Names: 4-1BB, ILA, CD137, TNFRSF9

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1509>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>