

PE/Cyanine5 anti-human CD137 (4-1BB)

Mouse IgG1, kappa
4B4-1



BioLegend
From Revvity

100 tests

Artikel-Nr.: B309808 | BioLegend | Hersteller-Nr.: 309808

253,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 4B4-1
Format: PE_CYANINE5
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 anti-human CD137 (4-1BB)

Clone: 4B4-1
Isotype: Mouse IgG1, kappa
Other_Names: 4-1BB, ILA, CD137, TNFRSF9
Size: 100 tests

Reactivity: Human
Format: PE_CYANINE5
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3909>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>