

Purified anti-human CD253 (TRAIL)

Mouse IgG1, kappa
RIK-2



BioLegend
From Revvity

100 µg

Artikel-Nr.: B308202 | **BioLegend** | Hersteller-Nr.: 308202

101,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: RIK-2
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD253 (TRAIL)

Clone: RIK-2
Isotype: Mouse IgG1, kappa
Other_Names: Apo-2 ligand, Apo-2L, TNFSF10, CD253
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2562>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>