

## Purified anti-human CD255 (TWEAK)

Mouse IgG3, kappa  
CARL-1



**BioLegend**  
From Revvity

100 µg

Artikel-Nr.: B308302 | BioLegend | Hersteller-Nr.: 308302

**156,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FC, IHC  
**Clone:** CARL-1  
**Format:** PURE  
**Isotype:** Mouse IgG3, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Purified anti-human CD255 (TWEAK)

Clone: CARL-1

Isotype: Mouse IgG3, kappa

Other\_Names: TNF-like weak inducer of apoptosis, Apo-3 ligand (Apo3L), Death receptor 3-L (DR3-L), TNFSF 12, CD255

Size: 100 µg

Reactivity: Human

Format: PURE

App\_Abbrev: FC, IHC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1333>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>