

PE anti-human CD255 (TWEAK)

Mouse IgG3, kappa
CARL-1



BioLegend
From Revvity

25 tests

Artikel-Nr.: B308305 | BioLegend | Hersteller-Nr.: 308305

108,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: CARL-1
Format: PE
Isotype: Mouse IgG3, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE anti-human CD255 (TWEAK)

Clone: CARL-1

Isotype: Mouse IgG3, kappa

Other_Names: TNF-like weak inducer of apoptosis, Apo-3 ligand (Apo3L), Death receptor 3-L (DR3-L), TNFSF 12, CD255

Size: 25 tests

Reactivity: Human

Format: PE

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1332>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>