

PE anti-human CD119 (IFN-gamma R alpha chain)

Mouse IgG1, kappa
GIR-208



BioLegend
From Revvity

100 tests

Artikel-Nr.: B308606 | BioLegend | Hersteller-Nr.: 308606

196,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: GIR-208
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE anti-human CD119 (IFN-gamma R alpha chain)

Clone: GIR-208

Isotype: Mouse IgG1, kappa

Other_Names: IFN-gammaR, IFN-gammaRalpha, IFN-gamma R1, CDw119

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2203>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>