

Purified anti-human CD210 (IL-10 R)

Rat IgG2a, kappa
3F9



BioLegend
From Revvity

100 µg

Artikel-Nr.: B308802 | **BioLegend** | Hersteller-Nr.: 308802

92,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 3F9
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 100 µg

Purified anti-human CD210 (IL-10 R)

Clone: 3F9

Isotype: Rat IgG2a, kappa

Other_Names: IL-10 Receptor, IL-10R

Size: 100 µg

Reactivity: Human, Cynomolgus, Rhesus

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1524>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>