

PE anti-human CD263 (DcR1, TRAIL-R3)

Mouse IgG1, kappa
DJR3



BioLegend
From Revvity

25 tests

Artikel-Nr.: B307005 | [BioLegend](#) | Hersteller-Nr.: 307005

113,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: DJR3
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE anti-human CD263 (DcR1, TRAIL-R3)

Clone: DJR3
Isotype: Mouse IgG1, kappa
Other_Names: TRAIL-R3, TRID, LIT, CD263, TNFRSF10C
Size: 25 tests

Reactivity: Human
Format: PE
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=637>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>