

Purified anti-human DR3 (TRAMP)

Mouse IgG1, kappa
JD3



BioLegend
From Revvity

100 µg
Artikel-Nr.: B307102 | **BioLegend** | Hersteller-Nr.: 307102

165,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: JD3
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human DR3 (TRAMP)

Clone: JD3

Isotype: Mouse IgG1, kappa

Other_Names: TNFRSF25, TRAMP, Apo-3, WSL-1, LARD, DDR3, TR3, TNFRSF12

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=786>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>