

PE anti-human DR3 (TRAMP)

Mouse IgG1, kappa
JD3



BioLegend
From Revvity

25 tests

Artikel-Nr.: B307105 | [BioLegend](#) | Hersteller-Nr.: 307105

113,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: JD3
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE anti-human DR3 (TRAMP)

Clone: JD3

Isotype: Mouse IgG1, kappa

Other_Names: TNFRSF25, TRAMP, Apo-3, WSL-1, LARD, DDR3, TR3, TNFRSF12

Size: 25 tests

Reactivity: Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=785>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>