

Purified anti-human CD261 (DR4, TRAIL-R1)

Mouse IgG1, kappa
DJR1



BioLegend
From Revvity

100 µg

Artikel-Nr.: B307202 | [BioLegend](#) | Hersteller-Nr.: 307202

165,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: DJR1
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD261 (DR4, TRAIL-R1)

Clone: DJR1
Isotype: Mouse IgG1, kappa
Other_Names: TRAIL-R1, Apo-2, CD261, TNFRSF10A
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=635>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>