

PE/Cyanine5 anti-human HLA-DR

Mouse IgG2a, kappa
L243



BioLegend
From Revvity

100 tests

Artikel-Nr.: B307608 | **BioLegend** | Hersteller-Nr.: 307608

235,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: L243
Format: PE_CYANINE5
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 anti-human HLA-DR

Clone: L243

Isotype: Mouse IgG2a, kappa

Other_Names: Major Histocompatibility Class II, MHC class II

Size: 100 tests

Reactivity: Human, Cynomolgus, Rhesus

Format: PE_CYANINE5

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=791>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>