

APC/Cyanine7 anti-human HLA-DR

Mouse IgG2a, kappa
L243



BioLegend
From Revvity

25 tests

Artikel-Nr.: B307617 | [BioLegend](#) | Hersteller-Nr.: 307617

115,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: L243
Format: APC_CYANINE7
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 25 tests

APC/Cyanine7 anti-human HLA-DR

Clone: L243

Isotype: Mouse IgG2a, kappa

Other_Names: Major Histocompatibility Class II, MHC class II

Size: 25 tests

Reactivity: Human, Cynomolgus, Rhesus

Format: APC_CYANINE7

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2863>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>