

Alexa Fluor® 700 anti-human HLA-DR

Mouse IgG2a, kappa
L243



BioLegend
From Revvity

25 µg
Artikel-Nr.: B307625 | **BioLegend** | Hersteller-Nr.: 307625

100,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: L243
Format: A700
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Alexa Fluor® 700 anti-human HLA-DR

Clone: L243
Isotype: Mouse IgG2a, kappa
Other_Names: Major Histocompatibility Class II, MHC class II
Size: 25 µg

Reactivity: Human, Cynomolgus, Rhesus
Format: A700
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3403>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>