

## Alexa Fluor® 700 anti-human HLA-DR

Mouse IgG2a, kappa  
L243



**BioLegend**  
From Revvity

100 µg

Artikel-Nr.: B307626 | BioLegend | Hersteller-Nr.: 307626

**212,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FC, ICFC  
**Clone:** L243  
**Format:** A700  
**Isotype:** Mouse IgG2a, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Alexa Fluor® 700 anti-human HLA-DR

Clone: L243

Isotype: Mouse IgG2a, kappa

Other\_Names: Major Histocompatibility Class II, MHC class II

Size: 100 µg

Reactivity: Human, Cynomolgus, Rhesus

Format: A700

App\_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3403>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>