

Alexa Fluor® 488 anti-human/mouse/rat PCNA

Mouse IgG2a, kappa
PC10



BioLegend
From Revvity

25 tests

Artikel-Nr.: B307909 | **BioLegend** | Hersteller-Nr.: 307909

100,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: PC10
Format: A488
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 25 tests

Alexa Fluor® 488 anti-human/mouse/rat PCNA

Clone: PC10

Isotype: Mouse IgG2a, kappa

Other_Names: Proliferating Cell Nuclear Antigen, DNA Polymerase delta Auxiliary Protein

Size: 25 tests

Reactivity: Human, Mouse, Rat

Format: A488

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3153>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>