

PE anti-human CD150 (SLAM)

Mouse IgG1, kappa
A12 (7D4)



BioLegend
From Revvity

25 tests

Artikel-Nr.: B306307 | **BioLegend** | Hersteller-Nr.: 306307

98,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: A12 (7D4)
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE anti-human CD150 (SLAM)

Clone: A12 (7D4)
Isotype: Mouse IgG1, kappa
Other_Names: SLAM, Signaling Lymphocyte Activation Molecule, IPO-3
Size: 25 tests

Reactivity: Human, Cynomolgus, Rhesus
Format: PE
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=590>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>