

Purified anti-human CD184 (CXCR4)

Mouse IgG2a, kappa
12G5



BioLegend
From Revvity

25 µg
Artikel-Nr.: B306501 | [BioLegend](#) | Hersteller-Nr.: 306501

67,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 12G5
Format: PURE
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-human CD184 (CXCR4)

Clone: 12G5
Isotype: Mouse IgG2a, kappa
Other_Names: CXCR4, Fusin
Size: 25 µg

Reactivity: Human, Cynomolgus, Rhesus
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=544>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>