

PE/Cyanine5 anti-human CD184 (CXCR4)

Mouse IgG2a, kappa
12G5



BioLegend
From Revvity

100 tests

Artikel-Nr.: B306508 | **BioLegend** | Hersteller-Nr.: 306508

252,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 12G5
Format: PE_CYANINE5
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 anti-human CD184 (CXCR4)

Clone: 12G5
Isotype: Mouse IgG2a, kappa
Other_Names: CXCR4, Fusin
Size: 100 tests

Reactivity: Human, Cynomolgus, Rhesus
Format: PE_CYANINE5
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=543>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>