

Purified anti-human CD235ab

Mouse IgG2b, kappa
HIR2 (GA-R2)



BioLegend
From Revvity

100 µg

Artikel-Nr.: B306602 | **BioLegend** | Hersteller-Nr.: 306602

82,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: HIR2 (GA-R2)
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD235ab

Clone: HIR2 (GA-R2)
Isotype: Mouse IgG2b, kappa
Other_Names: Glycophorin A/B, GPA/GPB, GYPA, GYPB
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=743>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>