

## PE anti-human CD235ab

Mouse IgG2b, kappa  
HIR2 (GA-R2)



**BioLegend**  
From Revvity

25 µg  
Artikel-Nr.: B306603 | **BioLegend** | Hersteller-Nr.: 306603

**83,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FC, ICFC  
**Clone:** HIR2 (GA-R2)  
**Format:** PE  
**Isotype:** Mouse IgG2b, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 25 µg

PE anti-human CD235ab

Clone: HIR2 (GA-R2)  
Isotype: Mouse IgG2b, kappa  
Other\_Names: Glycophorin A/B, GPA/GPB, GYPA, GYPB  
Size: 25 µg

Reactivity: Human  
Format: PE  
App\_Abbrev: FC, ICFC  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=741>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>