

PE/Cyanine5 anti-human CD235ab

Mouse IgG2b, kappa
HIR2 (GA-R2)



BioLegend
From Revvity

25 µg
Artikel-Nr.: B306605 | **BioLegend** | Hersteller-Nr.: 306605

105,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: HIR2 (GA-R2)
Format: PE_CYANINE5
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 25 µg

PE/Cyanine5 anti-human CD235ab

Clone: HIR2 (GA-R2)
Isotype: Mouse IgG2b, kappa
Other_Names: Glycophorin A/B, GPA/GPB, GYPA, GYPB
Size: 25 µg

Reactivity: Human
Format: PE_CYANINE5
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=742>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>