

Alexa Fluor® 488 anti-human TCR alpha/beta

Mouse IgG1, kappa
IP26



BioLegend
From Revvity

100 tests

Artikel-Nr.: B306712 | **BioLegend** | Hersteller-Nr.: 306712

212,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: IP26
Format: A488
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

Alexa Fluor® 488 anti-human TCR alpha/beta

Clone: IP26
Isotype: Mouse IgG1, kappa
Other_Names: alpha/beta TCR, TCR alpha/beta
Size: 100 tests

Reactivity: Human
Format: A488
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3373>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>