

PE/Cyanine5 anti-human CD106

Mouse IgG1, kappa
STA



BioLegend
From Revvity

100 tests

Artikel-Nr.: B305808 | [BioLegend](#) | Hersteller-Nr.: 305808

226,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: STA
Format: PE_CYANINE5
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 anti-human CD106

Clone: STA
Isotype: Mouse IgG1, kappa
Other_Names: VCAM-1, INCAM-110
Size: 100 tests

Reactivity: Human
Format: PE_CYANINE5
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=843>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>