

Biotin anti-human CD123

Mouse IgG1, kappa
6H6



BioLegend
From Revvity

100 µg
Artikel-Nr.: B306004 | **BioLegend** | Hersteller-Nr.: 306004

141,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 6H6
Format: BIOTIN
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Biotin anti-human CD123

Clone: 6H6
Isotype: Mouse IgG1, kappa
Other_Names: IL-3Ralpha, IL-3 Receptor alpha
Size: 100 µg

Reactivity: Human
Format: BIOTIN
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=576>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>