

PE/Cyanine7 anti-human CD123

Mouse IgG1, kappa
6H6



BioLegend
From Revvity

100 tests

Artikel-Nr.: B306010 | [BioLegend](#) | Hersteller-Nr.: 306010

249,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: 6H6
Format: PE_CYANINE7
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine7 anti-human CD123

Clone: 6H6
Isotype: Mouse IgG1, kappa
Other_Names: IL-3Ralpha, IL-3 Receptor alpha
Size: 100 tests

Reactivity: Human
Format: PE_CYANINE7
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2859>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>