

APC anti-human CD123

Mouse IgG1, kappa
6H6



BioLegend
From Revvity

25 tests

Artikel-Nr.: B306011 | **BioLegend** | Hersteller-Nr.: 306011

99,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 6H6
Format: APC
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

APC anti-human CD123

Clone: 6H6

Isotype: Mouse IgG1, kappa

Other_Names: IL-3Ralpha, IL-3 Receptor alpha

Size: 25 tests

Reactivity: Human

Format: APC

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3069>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>