

Purified anti-human CD147

Mouse IgG1, kappa
HIM6



BioLegend
From Revvity

100 µg

Artikel-Nr.: B306202 | [BioLegend](#) | Hersteller-Nr.: 306202

92,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: HIM6
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD147

Clone: HIM6
Isotype: Mouse IgG1, kappa
Other_Names: Neurothelin, Basigin, EMMPRIN
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=729>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>