

Alexa Fluor® 647 anti-human CD147

Mouse IgG1, kappa
HIM6



BioLegend
From Revvity

100 tests

Artikel-Nr.: B306210 | **BioLegend** | Hersteller-Nr.: 306210

212,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: HIM6
Format: A647
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

Alexa Fluor® 647 anti-human CD147

Clone: HIM6
Isotype: Mouse IgG1, kappa
Other_Names: Neurothelin, Basigin, EMMPRIN
Size: 100 tests

Reactivity: Human
Format: A647
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3368>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>