

Purified anti-human CD150 (SLAM)

Mouse IgG1, kappa
A12 (7D4)



BioLegend
From Revvity

100 µg
Artikel-Nr.: B306302 | **BioLegend** | Hersteller-Nr.: 306302

81,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: A12 (7D4)
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD150 (SLAM)

Clone: A12 (7D4)
Isotype: Mouse IgG1, kappa
Other_Names: SLAM, Signaling Lymphocyte Activation Molecule, IPO-3
Size: 100 µg

Reactivity: Human, Cynomolgus, Rhesus
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=591>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>