

Alexa Fluor® 647 anti-human CD80

Mouse IgG1, kappa
2D10



BioLegend
From Revvity

25 tests

Artikel-Nr.: B305215 | BioLegend | Hersteller-Nr.: 305215

95,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 2D10
Format: A647
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

Alexa Fluor® 647 anti-human CD80

Clone: 2D10
Isotype: Mouse IgG1, kappa
Other_Names: B7-1, B7, BB1
Size: 25 tests

Reactivity: Human
Format: A647
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3352>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>