

Purified anti-human CD83

Mouse IgG1, kappa
HB15e



BioLegend
From Revvity

25 µg
Artikel-Nr.: B305301 | **BioLegend** | Hersteller-Nr.: 305301

62,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: HB15e
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-human CD83

Clone: HB15e
Isotype: Mouse IgG1, kappa
Other_Names: HB15
Size: 25 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=683>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>