

Purified anti-human CD86

Mouse IgG2b, kappa
IT2.2



BioLegend
From Revvity

25 µg
Artikel-Nr.: B305401 | **BioLegend** | Hersteller-Nr.: 305401

62,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: IT2.2
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-human CD86

Clone: IT2.2
Isotype: Mouse IgG2b, kappa
Other_Names: B7-2, B70
Size: 25 µg

Reactivity: Human, Cynomolgus, Rhesus
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=780>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>