

PE/Cyanine5 anti-human CD86

Mouse IgG2b, kappa
IT2.2



BioLegend
From Revvity

100 tests

Artikel-Nr.: B305408 | [BioLegend](#) | Hersteller-Nr.: 305408

248,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: IT2.2
Format: PE_CYANINE5
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 anti-human CD86

Clone: IT2.2
Isotype: Mouse IgG2b, kappa
Other_Names: B7-2, B70
Size: 100 tests

Reactivity: Human, Cynomolgus, Rhesus
Format: PE_CYANINE5
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=779>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>