

Purified anti-Tubulin-gamma

Mouse IgG2b, kappa
14C11



BioLegend
From Revvity

25 µg
Artikel-Nr.: B629201 | **BioLegend** | Hersteller-Nr.: 629201

81,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ICC, WB
Clone: 14C11
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-Tubulin-gamma

Clone: 14C11
Isotype: Mouse IgG2b, kappa
Other_Names: gamma-tubulin
Size: 25 µg

Reactivity: Human, Mouse, Rat
Format: PURE
App_Abbrev: WB, ICC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3171>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>