

Purified anti-human TSLP

Mouse IgG1, kappa
13H8.3



BioLegend
From Revvity

500 µg

Artikel-Nr.: B512102 | BioLegend | Hersteller-Nr.: 512102

254,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA
Clone: 13H8.3
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-human TSLP

Clone: 13H8.3
Isotype: Mouse IgG1, kappa
Other_Names: Thymic stromal derived lymphopoietin
Size: 500 µg

Reactivity: Human
Format: PURE
App_Abbrev: ELISA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=4406>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>