

96er DeepWell Platte riplate, 1 ml, RW

U-Boden, Polypropylen
SBS Format



Ritter

Packung à 20 Stück

Artikel-Nr.: 742910 | Ritter | Hersteller-Nr.: 43001-0016

92,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Boden: U-Boden
Kategorie: 96er DeepWell Platte
Material: PP
Sterilität: unsteril
Verpackung: Packung à 20 Stück
Volumen: 1.0 ml
Well Form: rund

Riplate® RW 1.0 ml

- Erhöhte Wells verringern das Risiko der Kreuzkontamination und verbessern die Sicherheit beim Verschweißen oder Verschließen mit Klebefolie.
- Optimierte Kontur im SBS-Standard.
- Für Anwendungen im Thermocycler oder zum Auftauen gefrorener Proben im Wasserbad wurde die Riplate® plus konzipiert.
- Für Anwendungen mit extremen Temperaturbedingungen wie Tiefgefrieren oder Auftauen in der Mikrowelle.
- Der runde Boden garantiert eine hohe Probenrückgewinnung.
- Durch die chemische Resistenz ist Riplate® für sensible Diagnoseanwendungen geeignet.

Riplate® RW 1.0 ml

- *Raised wells reduce the risk of cross-contamination and increase the security of heat-sealing or sealing with adhesive films*

- *Optimized footprint according to the SBS standard*
- *For use in a thermocycler or for defrosting of samples in a hot-water bath, Riplate® plus was developed*
- *For applications with extreme temperature circumstances like freezing or defrosting in a microwave*
- *The U-shaped well bottom reduces the death volume*
- *Due to the high chemical resistance ideal for sensitive diagnostic applications.*