

96er DeepWell Platte riplate medio, 1.2 ml, RW

U-Boden, Polypropylen
SBS Format



Ritter
Packung à 32 Stück
Artikel-Nr.: 742917 | Ritter | Hersteller-Nr.: 43001-0101

154,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Boden: U-Boden
Eigenschaft: Medio
Kategorie: 96er DeepWell Platte
Material: PP
Verpackung: Packung à 32 Stück
Volumen: 0.5 ml
Well Form: rund

Riplate[®] medio

Bis zu 50% weniger Stapelhöhe

- Durch sich tangential berührende Wells benötigt Riplate[®] medio bis zu 50% geringere Stapelhöhe.
- Die erhöhten Ränder der einzelnen Wells lassen die Riplate[®] medio luftdicht verschweißen oder verkleben.
- Die alphanumerische Kennzeichnung ermöglicht die Registrierung und Identifizierung der einzelnen Proben.
- Neukonzipierte Wells ermöglichen eine gute Probendurchmischung.

Rplate[®] medio 1.0 ml

- Höhe: nur 27.0 mm
- Gesamtvolumen: 1.2 ml / Well
- Arbeitsvolumen: 50 – 1000 µl

Riplate[®] medio

- Due to the tangency of Riplate[®] medio, less than 50 % stacking size is needed.
- Higher well rims allow an airtight closure or sealing of Riplate[®] medio.
- The alphanumeric marking allows the registration and identification of single samples.
- New developed conical wells provide a good bioturbation.

Riplate[®] medio 1.0 ml

- Height: only 27.0 mm
- Total volume: 1.2 ml / well
- Working volume: 50 - 1000 µl / well