

96er DeepWell Platte riplate medio, 0.65 ml, RW

U-Boden, Polypropylen
SBS Format



Ritter

Packung à 40 Stück

Artikel-Nr.: 742918 | Ritter | Hersteller-Nr.: 43001-0102

199,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Boden: U-Boden

Eigenschaft: Medio

Kategorie: 96er DeepWell Platte

Material: PP

Verpackung: Packung à 40 Stück

Volumen: 0.5 ml

Well Form: rund

Riplate[®] medio

Bis zu 50% weniger Stapelhöhe

- Durch sich tangential berührende Wells benötigt Riplate[®] medio bis zu 50% geringere Stapelhöhe.
- Die erhöhten Ränder der einzelnen Wells lassen die Riplate[®] medio luftdicht verschweißen oder verkleben.
- Die alphanumerische Kennzeichnung ermöglicht die Registrierung und Identifizierung der einzelnen Proben.
- Neukonzipierte Wells ermöglichen eine gute Probendurchmischung.

Riplate[®] medio 0.5 ml

- Höhe: nur 20.5 mm
- Gesamtvolumen: 0.65 ml / Well
- Arbeitsvolumen: 20 - 500 µl / Well

Riplate® medio

- Due to the tangency of Riplate® medio, less than 50 % stacking size is needed.
- Higher well rims allow an airtight closure or sealing of Riplate® medio.
- The alphanumeric marking allows the registration and identification of single samples.
- New developed conical wells provide a good bioturbation.

Riplate® medio 0.5 ml

- Height: only 20.5 mm
- Total volume: 0.65 ml / well
- Working volume: 20 - 500 µl / well