

PE anti-human CD85d (ILT4)

Rat IgG2a, kappa
42D1



BioLegend
From Revvity

100 tests

Artikel-Nr.: B338706 | BioLegend | Hersteller-Nr.: 338706

254,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: 42D1
Format: PE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 100 Tests

PE anti-human CD85d (ILT4)

Clone: 42D1

Isotype: Rat IgG2a, kappa

Other_Names: CD85, ILT4, MIR10, Leukocyte Immunoglobulin-Like Receptor subfamily B member 2 (LILRB2), Leukocyte Immunoglobulin-like Receptor 2 (LIR-2)

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=5576>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>