

PerCP/Cyanine5.5 anti-human CD271 (NGFR)

Mouse IgG1, kappa
ME20.4



BioLegend
From Revvity

25 tests

Artikel-Nr.: B345111 | BioLegend | Hersteller-Nr.: 345111

141,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: ME20.4
Format: PERCP_CYANINE5.5
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PerCP/Cyanine5.5 anti-human CD271 (NGFR)

Clone: ME20.4

Isotype: Mouse IgG1, kappa

Other_Names: p75NTR, TNFRSF16, p75(NTR), Gp80-LNGFR, NGFR, Nerve Growth Factor Receptor

Size: 25 tests

Reactivity: Human

Format: PERCP_CYANINE5.5

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=7834>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>