

PE anti-mouse DcTRAIL-R1

Armenian Hamster IgG
mDcR1-3



BioLegend
From Revvity

100 µg

Artikel-Nr.: B133804 | BioLegend | Hersteller-Nr.: 133804

259,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: mDcR1-3
Format: PE
Isotype: Armenian Hamster IgG
Species of Origin: Armenian Hamster
Verpackung: 100 µg

PE anti-mouse DcTRAIL-R1

Clone: mDcR1-3

Isotype: Armenian Hamster IgG

Other_Names: TNFRH1, mSOB, Tumor necrosis factor receptor superfamily member 23, Tnfrsf23, Tumor necrosis factor receptor p60 homolog 1, TNF receptor family member SOB, Decoy TRAIL receptor 1, TNF receptor homolog 1

Size: 100 µg

Reactivity: Mouse

Format: PE

App_Abbrev: FC

Species of Origin: Armenian Hamster

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=5894>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>