

Purified anti-human RANTES (CCL5)

Mouse IgG2b, kappa
VL1



BioLegend
From Revvity

100 µg

Artikel-Nr.: B515502 | BioLegend | Hersteller-Nr.: 515502

123,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA, IHC, WB
Clone: VL1
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human RANTES (CCL5)

Clone: VL1

Isotype: Mouse IgG2b, kappa

Other_Names: CCL5, C-C motif chemokine 5, Small-inducible cytokine A5, T-cell-specific protein RANTES, SIS-delta, T cell-specific protein P228, TCP228, Eosinophil-chemotactic cytokine, EoCP

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: WB, IHC, ELISA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=6058>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>