

PE anti-mouse CD88 (C5aR)

Rat IgG2b, kappa
20/70



BioLegend
From Revvity

50 µg

Artikel-Nr.: B135805 | BioLegend | Hersteller-Nr.: 135805

117,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: 20/70
Format: PE
Isotype: Rat IgG2b, kappa
Species of Origin: Rat
Verpackung: 50 µg

PE anti-mouse CD88 (C5aR)

Clone: 20/70

Isotype: Rat IgG2b, kappa

Other_Names: C5a ligand, C5a anaphylatoxin receptor, C5AR, C5R1, C5AR1

Size: 50 µg

Reactivity: Mouse

Format: PE

App_Abbrev: FC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=6240>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>