

## Purified anti-human CD90 (Thy1)

Mouse IgG1, kappa  
5E10



**BioLegend**  
From Revvity

100 µg

Artikel-Nr.: B328102 | BioLegend | Hersteller-Nr.: 328102

**141,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** CyTOF®, FC, ICC, IHC-F, IP  
**Clone:** 5E10  
**Format:** PURE  
**Isotype:** Mouse IgG1, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Purified anti-human CD90 (Thy1)

Clone: 5E10  
Isotype: Mouse IgG1, kappa  
Other\_Names: Thy-1, Thy1  
Size: 100 µg

Reactivity: Human  
Format: PURE  
App\_Abbrev: FC, CyTOF®, IHC-F, ICC, IP  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:  
<https://www.biolegend.com/Default.aspx?ID=6664&productid=4110>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>