

Purified anti-human TSLP

Goat Polyclonal Ig
Poly5159



BioLegend
From Revvity

100 µg

Artikel-Nr.: B515901 | BioLegend | Hersteller-Nr.: 515901

280,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA Capture, WB
Clone: Poly5159
Format: PURE
Isotype: Goat Polyclonal Ig
Species of Origin: Goat
Verpackung: 100 µg

Purified anti-human TSLP

Clone: Poly5159
Isotype: Goat Polyclonal Ig
Other_Names: Thymic stromal lymphopoietin (TSLP)
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: ELISA Capture, WB
Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=8149>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>